

# TEST PAPERS

pH Test Paper-Rolls	38
pH Test Paper-Booklets	39
pH Test Paper-Strips	39
Litmus Paper-Booklets	40
Ion Test Papers	40
Chlorine Test Papers	40



## pH Test Paper-Booklets

- Narrow pH Unit Graduation: 0.2~0.4 (except UNIV)
- Can be used for severe pH measuring (except UNIV)
- Size: 7x70mm (20 leaves/booklet, 10 booklets/box)

Type		Catalog No.	Measuring Range pH
CR	Cresol red	07010010	0.4 ~ 2.0, 7.2 ~ 8.8
TB	Thymol blue	07010020	1.4 ~ 3.0, 8.0 ~ 9.6
BPB	Bromophenol blue	07010030	2.8 ~ 4.4
PB	Phenol blue	07010090	3.2 ~ 5.6
PP	Phenol purple	07010130	3.4 ~ 6.4
BCG	Bromocresol green	07010040	4.0 ~ 5.6
CPR	Chlorophenol red	07010100	5.0 ~ 6.6
MR	Methyl red	07010050	5.4 ~ 7.0
BCP	Bromocresol purple	07010140	5.6 ~ 7.2
BTB	Bromothymol blue	07010060	6.2 ~ 7.8
PR	Phenol red	07010150	0.0 ~ 1.6, 6.6 ~ 8.2
AZY	Alizarin yellow	07010070	10.0 ~ 12.0
ALB	Alkali blue	07010080	11.0 ~ 13.6
UNIV	Universal	07010120	1.0 ~ 11.0
No.20	MR & BTB	07010110	5.0 ~ 8.0



## pH Test Paper-Strips

- Packed in Plastic Bottle
- Size: 7mm x 40mm  
(PB: 200 strips/bottle, BTB and AZY: 300 strips/bottle)

Type		Catalog No.	Measuring Range pH
PB	Phenol blue	08001040	3.2 ~ 5.6
BTB	Bromothymol blue	08001100	6.2 ~ 7.8
AZY	Alizarin yellow	08001120	10.0 ~ 12.0



## Litmus Paper-Booklets

- Detecting Acid (Blue) and Alkali (Red)
- High sensitivity

Type	Catalog No.	Size (mm)	Package
Litmus Paper (Red)	07020020	7 x 70	20 leaves/booklet, 10 booklets / box
Litmus Paper (Blue)	07020010	7 x 70	20 leaves/booklet, 10 booklets / box



## Ion Test Papers

- Detecting Ion and its concentration
- Measuring pH

Type	Measuring for	Measuring Range	Color	pH	Catalog No.	Package (strips/box)
Almi Check	Al <sup>3+</sup>	0 ~ 100mg/L	Peach to Red	3 ~ 6	07030010	100
Copper Check	Cu <sup>2+</sup>	0 ~ 500mg/L	Bark to Black and Blue	3 ~ 5	07030020	100
Copper Check A	Cu <sup>+</sup> , Cu <sup>2+</sup>	0 ~ 50mg/L	White to Purple	1 ~ 6	07030030	100
Chrome Check	Cr <sup>6+</sup>	0 ~ 200mg/L	Light peach to Purple	1 ~ 4	07030041	25
Chrome Check A	Cr <sup>6+</sup>	0 ~ 50mg/L	Light peach to Purple	1 ~ 9	07030051	25
Iron Check	Fe <sup>2+</sup>	0 ~ 1,000mg/L	White to Red	1 ~ 6	07030060	100
Nickel Check	Ni <sup>2+</sup>	0 ~ 1,000mg/L	Grey to Red	2 ~ 7	07030070	100
Silver Check	Ag <sup>+</sup> in fixing Solution	0 ~ 10g/L	Yellow to Dark Brown	-	07030080	100
Cyanogen Check	CN <sup>-</sup>	0 ~ 1,000mg/L	White to Bark	7 ~ 9	07030100	100

## Chlorine Test Papers

- Detecting the concentration of residual chlorine in sodium hypochlorite solution (NaClO) or hypochlorous acid (HClO)
- Two types are available

Type	Catalog No.	Size (mm)	Package (strips/box)
Chlorine Test Paper (10 ~ 50ppm)	08010010	7 x 40	300
Chlorine Test Paper (25 ~ 200ppm)	08010020	7 x 40	300

Note: Cannot be used for the detection of Chlorine Ion.

