

## 1.5ml Tube

THE OPTICAL CLARITY AND THE STRENGTH OF THE HIGH IMPACT POLYSTYRENE MAKES THIS TUBE IDEAL FOR SEROLOGY AND HAEMATOLOGY USE.

### POLYPROPYLENE TUBE WITH POLYETHYLENE CAP

**P4411** TUBE ONLY

**P4411UU** UNLABELLED, NATURAL PLUG CAP SEPARATE

### POLYSTYRENE TUBE WITH POLYETHYLENE CAP

**S4411** TUBE ONLY

**S4411UU** UNLABELLED, NATURAL PLUG CAP SEPARATE

**DIMENSIONS**  
44mm HIGH  
11mm WIDE  
(ACTUAL SIZE  
EXTERNAL)



INNER PACK 1000 TO BAG  
OUTER CARTON QUANTITY 5000  
CARTON WEIGHT 9.0kg  
CARTON DIMENSIONS 470 x 310 x 475mm

## 5ml Tube

### POLYPROPYLENE TUBE WITH POLYETHYLENE CAP

**P7611** TUBE ONLY

**P7611UU** UNLABELLED, NATURAL PLUG CAP SEPARATE

### POLYSTYRENE TUBE WITH POLYETHYLENE CAP

**S7611** TUBE ONLY

**S7611UU** UNLABELLED, NATURAL PLUG CAP SEPARATE

**DIMENSIONS**  
76mm HIGH  
11.5mm WIDE  
(ACTUAL SIZE  
EXTERNAL)



INNER PACK 500 TO BAG  
CARTON QUANTITY 3000  
CARTON WEIGHT 10.5kg  
CARTON DIMENSIONS 470 x 310 x 475mm

## Cap to suit Above

**E7611C-02** PLUG STYLE

