

# IKA®

## RCT Basic Promotion

# INTEGRATED SAFETY

Experience safety first hand with the RCT basic from IKA®

### Features

- Integrated temperature control
- Supplied with PT1000 temperature sensor
- Exact temperature and speed setting via digital display, even when switched off
- Digital display of set safety temp. limit
- Hot top indicator - hot surface warning to prevent burns!
- Digital error code display
- Adjustable safety circuit of heating plate temperature (50 to 360 degrees C)
- Highly polished aluminium heating plate for optimal heat transfer
- Supplied with protective cover
- Bushing according to DIN 12878 for connecting a contact thermometer (not supplied with delivery)

**\$599.00** ex-gst

Product Code: 38100001



Designed  
to work perfectly

Technical Data	
Max. Stirring Quantity (H2O)	20L
Speed Display	Digital
Speed Range	50 to 1500rpm
Max. Magnetic Bar	80 x 10mm
Temperature Range	RT to 310 degrees C
Heating Plate	Aluminium 135mm D


For the full technical data, please contact Labtek.

All pricing excludes GST. Specials expire: 31/12/11

Ph: **1300 881 318**  
Fx: **1300 881 513**

PO Box 5316 Brendale Q 4500  
Email: [sales@labtek.com.au](mailto:sales@labtek.com.au)  
Website: [www.labtek.com.au](http://www.labtek.com.au)

**Brisbane • Sydney • Melbourne • Townsville**



Labtek Laboratory & Scientific Supplies

Photographs within brochure are for illustrative use only