

ADAPTERS - "distillation"

➤ Made from borosilicate glass 3.3 according to DIN 12594.

➤ Suitable for distillation assemblies for connecting flasks to condenser.

➤ Upper socket designed to accept thermometer either by screw thread adapter or by thermometer socket.

➤ With side cone at 75° angled to the center socket.

**ADAPTERS - "distillation" - "claisen head"**

➤ Made from borosilicate glass 3.3 according to DIN 12594.

➤ Suitable for distillation assemblies for connecting flasks to condenser.

➤ Upper socket designed to accept thermometer either by screw thread adapter or by thermometer socket.

➤ With side cone at 75° angled to the center socket.

➤ Claisen head offers one more extra socket.



catalogue number	socket NS size	cone NS size	angled cone NS size	pack quantity
065.45.001	NS 14/23	NS 14/23	NS 14/23	1 piece
065.45.002	NS 14/23	NS 19/26	NS 19/26	1 piece
065.45.005	NS 14/23	NS 29/32	NS 29/32	1 piece

catalogue number	socket NS size	cone NS size	angled cone NS size	pack quantity
065.46.001	NS 14/23	NS 14/23	NS 14/23	1 piece
065.46.002	NS 14/23	NS 19/26	NS 19/26	1 piece
065.46.005	NS 14/23	NS 29/32	NS 14/23	1 piece
065.46.006	NS 14/23	NS 29/32	NS 29/32	1 piece

ADAPTERS - "delivery"

➤ Made from borosilicate glass 3.3. Suitable for distillation assemblies with vacuum connections.

➤ Available 2 types: ➤ plain, bend with angle 105°
➤ with vent, bend with angle 105°

adapters plain	adapters with vent	socket NS size	cone NS size	pack quantity
065.32.014	065.33.014	NS 14/23	NS 14/23	1 piece
065.32.029	065.33.029	NS 29/32	NS 29/32	1 piece

**ADAPTERS - "receiver"**

➤ Made from borosilicate glass 3.3. Suitable for distillation assemblies with vacuum connections.

➤ Available 2 types: ➤ straight
➤ bend with angle 105°

adapters straight	adapters bend	socket NS size	cone NS size	pack quantity
065.39.014	065.40.014	NS 14/23	NS 14/23	1 piece
065.39.029	065.40.029	NS 29/32	NS 29/32	1 piece

