



ADAPTERS - "splash head"

- Made from borosilicate glass 3.3. Suitable to be used in distillation assemblies to prevent raw liquid step over from flask to condenser.
- Available 2 types:
 - straight
 - bend with angle 75°

straight adapters	bend adapters	socket NS size	cone NS size	pack quantity
065.47.001	065.48.001	NS 14/23	NS 14/23	1 piece
065.47.003	065.48.005	NS 29/32	NS 29/32	1 piece

DRYING TUBES

- Made from borosilicate glass 3.3. In compliance with DIN 12616 standards.
- Available 2 types:
 - bend with angle 90°
 - U form

bend	U form	U form with stopcock	pack quantity
065.60.100	065.62.100	065.63.100	1 piece
065.60.150	065.62.130	065.63.130	1 piece
	065.62.180	065.63.180	1 piece



WASH BOTTLES - "gas"

- Made from borosilicate glass 3.3.
- Thick base with excellent stability.
- Smooth NS 29/32 neck are in compliance with DIN 12596 standards.



WASH BOTTLE HEADS

- Made from borosilicate glass 3.3.
- Smooth NS 29/32 neck are in compliance with DIN 12242 standards.
- Supplied with or without sintered glass support (of porosity 1).
- Sintered glass support which enables efficient gas dispensing.



catalogue number	bottle volume	socket NS size	pack quantity
065.66.125	125 ml	NS 29/32	1 piece
065.66.250	250 ml	NS 29/32	1 piece
065.66.500	500 ml	NS 29/32	1 piece
065.66.901	1000 ml	NS 29/32	1 piece

catalogue number	product description	pack quantity
065.67.001	without frit disk - with glass side arm	1 piece
065.67.002	without frit disk - with P.P side arm	1 piece
065.67.011	with frit disk - with glass side arm	1 piece
065.67.022	with frit disk - with P.P side arm	1 piece